Your practice is *our* inspiration.™

SonicFill™ Composite

Sonic-Activated Nano-Hybrid Composite

Composite Material

Instructions For Use

SonicFill

Nanohybrid Composite Restorative

SonicFill is a light-cured, low-shrink, resin-based, dental restorative designed for direct placement. It is indicated for all cavity classes in posterior teeth. SonicFill is used in combination with the SonicFill Handpiece, and offers sonically activated delivery. This delivery produces a significant drop in viscosity upon extrusion allowing intimate adaptation of the composite to the cavity walls. When the cavity is filled and the Handpiece is deactivated, the viscosity of the restorative resin increases retaining desired physical properties. The non-sticky, non-slumping nature of the material allows for quick, easy shaping and sculpting. These benefits, along with a high depth of cure, allow a technique in which a cavity up to 5mm in depth can be filled and cured in a single bulk increment.

Storage Conditions and Shelf Life:

SonicFill should be stored at ambient temperature. The expiration date is marked on packaging. Do not use after expiration date. Every effort should be made to protect the composite from visible light prior to use.

PLEASE NOTE THE FOLLOWING:

The Unidose® tips are designed for SINGLE PATIENT USE ONLY. Do not re-cap and/or re-use the Unidose® tip once material has been dispensed for that patient. This is to prevent cross-contamination between patients.

PRIOR TO PLACEMENT --

RECOMMENDATIONS ON PROPER BONDING

- Isolation throughout adhesive steps and composite placement is important. Rubber dam is ideal.
- Please closely follow bonding agent directions for use.
- Please take care to ensure that your air line is free of oil and other contaminants.

SonicFill Handpiece

The SonicFill Handpiece is packaged with a Use and Care guide. To receive the best use from your Handpiece, read and follow the use and care instructions included with the device.

PLACEMENT OF

SonicFill

1. Select the desired shade. If Class II, place matrix of choice.
2. Verify air pressure of dental unit is at least 36 PSI or 2,5 bar. If the unit is set to a higher pressure, it will be automatically regulated to 36 PSI or 2,5 bar inside the instrument.
3. Accurately position the Handpiece on the MULTIflex coupling and press it firmly until securely (audibly) locked.
4. Remove red protective covers from composite chamber and dispensing tip of Unidose capsule.
5. Insert Unidose capsule and properly attach to Handpiece by screwing it on clockwise with finger pressure.
6. Dispense with a Disk-type (on /off) foot control. The dispensing rate/speed is set with the switch at the bottom of the Handpiece. Setting 5 is the fastest speed; setting 1 is the slowest. Set desired speed at base of Handpiece. For your first several uses, you may wish to place the speed at setting "3" until you become familiar with the flow rate. Afterwards, the different dispensing rate/speed setting should be used based on indications and clinical situation.
Dispensing with an infinitely variable foot control (e.g. KaVo multifunctional foot control): Infinitely variable foot controls are able to set the dispensing rate/speed via the foot pedal position. For example, the KaVo foot control will deliver maximum volume with a pedal position to the right, where as the position to the left will reduce the delivered volume. In this case set the switch at the bottom of the Handpiece to "5".
7. Place Unidose® tip at the deepest portion of the preparation to avoid trapping air.
8. Activate SonicFill Handpiece by depressing foot pedal and fill entire cavity (up to 5 mm). It is recommended to use full air pressure when activating the handpiece. Then use the dispensing rate switch or the foot control to adjust the flow rate and speed of the composite. In a deep cavity (>5mm) or when placing core buildup in a deep pulp chamber, 2 separately cured increments are recommended.
9. After placement, define anatomy using hand instrument(s).
10. Light Cure *.

11. Adjust occlusion, finish and polish in the usual manner.
12. Remove tip by unscrewing it counterclockwise with finger pressure.
13. Between patients, follow the infection control and maintenance guidelines specified in the SonicFill Handpiece Use and Care Guide.

Recommended Cure Times: Demi/Demi Plus - 20 seconds, L.E.Demetron II - 20 seconds; or Optilux 501 in Boost mode - 20 seconds, Ramp Mode - 40 seconds, or Regular Mode - 40 seconds. For all other lights, see manufacturer's recommendation.

- * In the posterior, light cure the recommended time from the occlusal, remove the matrix and cure again from the buccal and lingual.
In a Class I, additional cure is still recommended from the facial and lingual.

CAUTION: Uncured methacrylate resin may cause contact dermatitis and damage the pulp. Avoid contact with skin, eyes and soft tissue. Wash thoroughly with water after contact.

Limited Warranty - Limitation of Kerr's Liability

Kerr's technical advice, whether verbal or in writing, is designed to assist dentists in using Kerr's product. The dentist assumes all risk and liability for damages arising out of the improper use of Kerr's product. In the event of a defect in material or workmanship, Kerr's liability is limited, at Kerr's option, to replacement of the defective product or part thereof, or reimbursement of the actual cost of the defective product. In order to take advantage of this limited warranty, the defective product must be returned to Kerr. In no event shall Kerr be liable for any indirect, incidental, or consequential damages.

EXCEPT AS EXPRESSLY PROVIDED ABOVE, THERE ARE NO WARRANTIES, BY KERR, EXPRESS OR IMPLIED, INCLUDING WARRANTIES WITH RESPECT TO DESCRIPTION, QUALITY, OR FITNESS FOR A PARTICULAR PURPOSE.

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Kerr Corporation
Kerr Italia, S.r.l.