

Demi only



921677

2mm

Curved

For layer curing of posterior restorations, inside gingival box preparation, and anterior tacking of laminate veneers.

Demi only



921676

4mm

Turbo Curved

For preventing overlapping on the gingiva in Class V.



20945

2mm

Curved

For layer curing of posterior restorations, inside gingival box preparation, and anterior tacking of laminate veneers.



21185

4mm

Turbo Curved

For preventing overlapping on the gingiva in Class V.



921535

8mm

Curved Amber

For general and posterior curing.



21020

8mm

Curved Turbo

Increases efficiency of light output and reduces overall curing time.



20941

8mm

Curved

For improved access in general and posterior curing.



952213

8mm

Turbo+

Increases efficiency of light output and reduces overall curing time.



20898

11mm

Curved

For improved access in general and posterior curing.



921551

11mm

Curved Turbo

Increases efficiency of light output and reduces overall curing time.



20812

13mm

Curved

For treating entire occlusal surface of molars.

Demetron Accessories



20879

Light Guide Holder

Holds six accessory light guides.



21042

Optics Maintenance Kit

After numerous autoclave cycles, brown discoloration may accumulate on polished ends of the light guide. Avoid by polishing ends every 50 cycles. Do not use to remove composite.

Your practice is *our* inspiration.™



Demetron

Light Guides

kerrdental.com 800.KERR.123

0309 5k